

SNAP *Giardia* Test

Test accuracy

Sensitivity and specificity of the SNAP® *Giardia* Test


Analyte	Reference method result*	Device result		Total	Sensitivity (95% CL)
		+	–		Specificity (95% CL)
<i>Giardia</i>	+	74	4	78	95% (87%–98%)
	–	1	144	145	99% (96%–100%)

CL = confidence limit

*Reference method was direct immunofluorescence microscopy.

Detecting *Giardia* infection

A comparison of in-house methods vs. reference laboratory methods



In a study, 617 fecal samples were tested for *Giardia* in the normal course of clinical practice.¹ Results from fecal flotation or fecal smears were recorded, and the samples were then forwarded for further reference testing. The results indicated poor performance of current in-house microscopy testing compared to enzyme-linked immunosorbent assay (ELISA) and direct immunofluorescence microscopy (DIFM) reference methods.

The SNAP *Giardia* Test and in-clinic microscopy test results are shown as percentages of samples that were identified as *Giardia*-positive by laboratory ELISA and DIFM methods.

The SNAP *Giardia* Test is finding positives that microscopy misses.¹

Reference

1. Groat R, Monn M, Flynn L, Curato J. Survey of clinic practices and testing for diagnosis of *Giardia* infections in dogs and cats. Presented at: 2003 ACVIM Forum; June 4–8; Charlotte, NC.

For detailed sensitivity and specificity information, please reference the product insert.

© 2021 IDEXX Laboratories, Inc. All rights reserved. • 09-80938-02

All ®/TM marks are owned by IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

The IDEXX Privacy Policy is available at idexx.com.

