Hearworm Antigen-Anaplasma-Leishmania infantum-Ehrlichia Antibody Test Kit

For veterinary use only.

English version 

SNAP* Leish 4Dx* Test

In vitro diagnostic for the detection of *Dirofilaria immitis* antigen, antibody to *Anaplasma phagocytophylum*, antibody to *Anaplasma platys*, antibody to *Leishmania infantum*, antibody to *Ehrlichia canis*, and antibody to *Ehrlichia ewingii* in canine serum, plasma, or whole blood.

Precautions and warnings

- All wastes should be properly decontaminated prior to disposal.
- Do not mix components from kits with different serial numbers.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- Male-only heartworm infections typically produce antigen levels that are below the detection capability of this antigen test.
- Refer to country specific Safety Data Sheet for regional hazard identification.

Storage

- Store at 2–8°C until expiration date.
- SNAP devices and reagents can be stored at room temperature (18–25°C) for 6 months or until the expiration date, whichever occurs first.
- After SNAP devices and reagents are removed from 2–8°C for more than 24 hours, the expiration date is 6 months or the printed expiration date, whichever occurs first. If the 6 month expiration date occurs prior to the printed expiration date, record the new date in the space provided on the kit.

Kit components

Item	Reagents	Quantity
1	1 bouteille Anti-D. <i>immitis</i> /A. <i>phagocytophylum</i> /A. <i>platys</i> /L. <i>infantum</i> /E. <i>canis</i> /E. <i>ewingii</i> : HRPO conjugué (conservé avec ProClin™ 150)	7,0 mL
2	SNAP dispositif	5, 15, 30
Reagents contained in each device:		
Wash solution (conservé avec ProClin™ 150)	0,4 mL	
Substrate solution	0,6 mL	
Other components: transfert pipettes, sample tubes, and reagent rack		

Conjugat - H317/H412/P261/P273/P280/P330+P310/P362+P364/P501: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapours. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose of contents/container in accordance with local/regional/international regulations.

Wash Solution - H317/H412/P261/P273/P280/P330+P310/P362+P364/P501: May cause an allergic skin reaction. Causes serious eye irritation. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapours. Avoid release to the environment. Wear eye protection/face protection/protection glasses. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

Sample information

- Samples must be at room temperature (18–25°C) before beginning the test procedure.
- Serum, plasma, or anti-coagulated whole blood (e.g., EDTA, heparin), either fresh or stored at 2–8°C for up to one week, can be used.
- For longer storage, serum, or plasma can be frozen (-20°C or colder) and then reconstituted before use.
- Hemolyzed or lipemic samples will not affect test results.

Test procedure

1. Stored in a refrigerator, allow all components to equilibrate at room temperature (18–25°C) for 30 minutes before use. **Do not heat.**

2. Using the pipette provided, dispense 3 drops of sample into a new sample tube.

3. Holding the bottle vertical, add 4 drops of conjugate to the sample tube.

4. Cap the sample tube and mix it thoroughly by inverting it 3–5 times.

5. Place the device on a horizontal surface. Add the entire contents of the sample tube to the sample well, being careful not to splash the contents outside of the sample well.

- The sample will flow across the result window, reaching the activation circle in 30–60 seconds. Some sample may remain in the sample well.

6. When color EBSI appears in the activation circle, push the activator firmly until it is flush with the device body.

- Note: Some samples may not flow to the activation circle within 60 seconds, and therefore, the circle may not turn color. In this case, press the activator after the sample has flowed across the result window.

7. Read test result at 8 minutes.

- Note: The positive control may develop sooner, but results are not complete until 8 minutes.

Interpreting the test results

Positive result

Any color development in the sample spots indicates the presence of *Dirofilaria immitis* antigen, A. *phagocytophylum* antibody, A. *platys* antibody, L. *infantum* antibody, E. *canis* antibody, or E. *ewingii* antibody in the sample.

Positive control → A. *phagocytophylum*/A. *platys*
E. *canis*/E. *ewigii* → A. *phagocytophylum*/A. *platys*
L. *infantum* → D. *immitis*

Notes:

- The A. *phagocytophylum*/A. *platys* spot cannot differentiate between the two species: a positive result indicates presence of antibodies to A. *phagocytophylum* and/or A. *platys*.
- The E. *canis*/E. *ewigii* spot cannot differentiate between the two species: a positive result indicates presence of antibodies to E. *canis* and/or E. *ewigii*.
- In a low percentage of samples (0.027% as reported), interfering substances in the patient's blood can cause all spots on the device to react as positive. In this event, the sample should be re-assayed as plasma or serum to reduce the likelihood of interference.

Negative result

- Only positive control spot develops color.

Invalid results

- Background—if the sample is allowed to flow past the activation circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No color development—if the positive control does not develop color, repeat the test.

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Conjugat - H317/H412/P261/P273/P280/P330+P313/P362+P364/P501: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapours. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose of contents/container in accordance with local/regional/international regulations.

Wash Solution - H317/H412/P261/P273/P280/P330+P313/P362+P364/P501: May cause an allergic skin reaction. Causes serious eye irritation. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapours. Avoid release to the environment. Wear eye protection/face protection/protection glasses. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

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